

Meet the Family

Safety & Efficiency

XP Focus

Nd:YAG 10w

- Onychomycosis
- Vacular lesions
- Acne treatment
- Anti aging
- Sagging skin improvements
- Hair removal
- Warts



XP-2 Focus

Nd:YAG 20w

- Endovascular
- Hyperhidrosis
- Laser Lipo
- Varicose veins
- Vacular lesions
- Onychomycosis
- Hair removal
- Acute acne
- Skin tightening

SP Spectro

Nd:YAG 35w & Er(3J)

- Ultimate safety
- Erb applications
- Full functional scanning
- Skin resurfacing
- Novel Turbo & V-smooth
- Complete anti-aging
- Fractional Er & Nd:YAG
- Endovascular
- Laser Lipo



SP Dynamis

Nd:YAG 80w & Er: YAG 3J

- Very powerful
- Ablative non Co² heat
- Skin resurfacing
- Skin tightening (non ablative)
- Fractional Er & Nd:YAG
- Hair removal (large scanner)
- Onychomycosis
- Endovascular
- Hyperhidrosis
- Laser Lipo

The Highest Performance - Best Made Laser Systems in the world

Meet the Family

Safety & Efficiency

XP Max

Nd:YAG 130w

- High Speed HR
- Scanner 42 cm²
- Largest spot size
- FRAC3[®] skin rejuvenation
- Wrinkle reduction
- Wart treatment
- Acne treatment
- Advanced VSP Technology (Square pulse)



QX MAX

Nd:YAG Q-Switch

- Q-switch laser
- Uniform beam profile
- Five laser source in one system
- Treats all pigmented and tattoo colours
- FRAC3[®] rejuvenation treatment
- Hair removal, acne
- No Cr 3 ions in rods (very stable)



Lightwalker ST-E

Er: 8w

- Dental Laser Er: 8w
- Variable pulse
- Square pulse
- Laser Endo (PIPS[®])ready



Lightwalker DT

Nd:YAG 8w & Er: 8w

- Tipless hand piece
- Er: Touch white
- Compressed air onboard
- Endodontics
- Whitening with Er
- PIPS[®] ready

Technolase cc

Tel: 012 349 1750 • Fax: 012 349 1752

Email: technola@mweb.co.za

www.technolase.co.za